
DentoNet
DAILY ORAL CARE

**Designing toothbrushes
from a timeless perspective.**

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
Spazzolificio Piave Spa started its activity over 150 years ago manufacturing all kinds of brushes (from hairbrushes to toothbrushes) and its main business was contract manufacturing for the Italian Army.

In 1970 the company was taken over by the Zecchin family, who transferred the activity to its present location, Onara. Here new models and brands were created for the distribution on national and international markets.

Research and design in the 80's improved the Piave toothbrushes developing more functional products with top quality raw material, creating many iconic products for the Italian and other prestigious markets.

In the 90's, with the launch of Silver Care, the antibacterial concept was introduced for the first time ever for this type of product. This innovative toothbrush was made possible due to a simple but effective solution: the naturally antibacterial silver coated replaceable head.

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Designing toothbrushes from a timeless perspective.

We offer oral hygiene solutions, ensuring efficient and functional products, developed in close cooperation with dental experts. As a constant process, we will continue to satisfy our customers without neglecting our planet. Our products are designed and manufactured in Italy.

All production is carried out in our factory in Onara ITALY, a unique feature for a company manufacturing toothbrushes and oral care accessories.

Made in Italy is prestigious. These words make every Italian manufacturer proud. It is based on creativity matched with quality and design. Our creativity and skills are properly valued on international and local markets.

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The Piave toothbrush company uses a solar panel system that independently produces clean and inexhaustible energy, avoiding the use of fossil fuels and consequently air pollution. We are actively involved in the development of new eco-friendly technologies, which is why we use eco sustainable materials wherever possible.

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In order to offer up-to-date models, in terms of quality and design, PIAVE uses the most advanced technology and closely cooperates with top oral care experts (Italian dentists and the Dentistry School of the University in Milano).

As further confirmation of the high quality standards, the company, is certified SO 9001-2015. The quality standards and Product Sterilization during the packaging process, make for top quality products.

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
Thanks to numerous clinical tests, a team of researchers confirms that for each pathology there is the right toothbrush, as well as for any dental arch conformation.

THE HEADS

The anatomical and rounded head, thinner at the tip, allows you to easily reach the less accessible parts of the mouth without causing trauma to the gums. The medium bristles with a rounded tip are obtained from renewable and natural sources.

THE HANDLE

The handle helps to correctly position the toothbrush and facilitate the essential rotary movements to implement the recommended brushing technique.

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
SOFT BRISTLES for gentle plaque removal and effective gum protection.

MEDIUM BRISTLES for more decisive action on bacterial plaque.

HARD BRISTLES ensure an anti-tartar action.

The dome-cut of the bristles allows the higher central rows to penetrate gently under the gum line. The shorter outer bristles effectively clean the dental surface. The perfect end-rounding of the tips avoids damage to tooth enamel and gums.

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HOW TO BRUSH


- Use a small amount of toothpaste.
- Brush the outer surface of all teeth using circular movements.
- Brush the inner tooth and chewing surfaces with perpendicular movements.
- Brush teeth for at least 2 minutes.
- Rinse with water.

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THE CORRECT BRUSHING TECHNIQUE

Place the toothbrush at 45 degrees, half on the gum, half on the tooth, and clean with small circular movements, from the gum towards the tooth and not vice-versa (modified Bass technique). Replace your toothbrush every 3 months.

Art. 4320
Morbido/Soft
EAN: 8009315043209
Art. 4321
Medio/Medium
EAN: 8009315043216
Art. 4322
Duro/Hard
EAN: 8009315043223

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This toothbrush is specific for those suffering from **dental sensitivity**.

The **soft microfine bristles** in fact, enter more easily within the interdental spaces and along the gum line, ensuring a more effective cleaning.

The **extreme softness** of the tips is ideal in cases of delicate gums, sensitive teeth or paradontal pathologies.

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HEALTHY GUMS, HEALTHY MOUTH

Gums protect and support teeth and tissues that attach them to the jaw. When your gums are not healthy there is a danger of losing teeth and altering your overall health.

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V-cut bristles specifically for effective cleaning of **orthodontic braces**.

The soft bristles, guarantee plaque removal without damaging tooth enamel.

The **compact head** allows for easy access to all points of the dental arch.

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ORAL HYGIENE WITH ORTHODONTIC BRACES

It should preferably be carried out scrupulously and constantly. Use a specific toothbrush with which it is possible to reach even the most hidden positions. In this case, brushing your teeth after every meal is even more important because bacterial plaque can build up more easily.

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
The travel toothbrush is ideal, for continued oral hygiene even when away from home. It may be dismantled thus attaining a compact structure that protects the bristles and is guaranteed to be as ergonomic as a traditional toothbrush.

Toothpaste with advanced formula: 5-action formula which fights cavities, tartar and gum inflammation. The strong mint flavour helps prevent halitosis. For adult use only.

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INSTRUCTIONS FOR USE

1. Wind about 40 cm of floss around your fingers and pull until taut.
2. Gently run the floss through the interdental spaces.

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Dental floss removes plaque that adheres to teeth and gums. Polishes tooth surfaces and prevents bad breath. Daily flossing decreases the possibility of having periodontal disease or tooth decay.

Dental floss with **fluoride** allows for the removal of bacterial plaque, and at the same time strengthens the tooth enamel, leaving a pleasant **mint**

flavour. The **extra-flat, waxed floss** allows for the removal of bacterial plaque in very narrow interdental spaces, assuring that it slides easily between teeth


Art. 102
Filo Fuoro + Menta
EAN: 8009315001025
Art. 108
Filo Extra Piatto
EAN: 8009315001087

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These **disposable interdental brushes** are gentle and easy to use. The soft rubber bristles, slide comfortably between teeth to remove plaque and food particles. They are ideal for cleaning between teeth and a great choice for very narrow interdental spaces, bridges, orthodontic braces or delicate gums.

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INTERDENTAL CLEANING

Using interdental brushes helps to prevent gum disease. An interdental brush removes plaque and bacteria that build up between teeth, that can irritate gums and cause swelling. Dentists recommend regular inter-dental cleaning.

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
DISPOSABLE FLOSSER WITH INTEGRATED DENTAL FLOSS EXTRA RESISTANT

Thanks to the practical and comfortable grip, maneuvering the flosser between the teeth becomes simple: ready to use, it effectively removes bacterial plaque in the interdental spaces and under the gum line. It is a perfect substitute for dental floss when away from home. The resealable bag maintains the integrity of the product.

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
Interdental brushes are the ideal instruments for effective cleaning of interdental spaces, orthodontic braces, crowns and bridges. They are easy to use and, thanks to their comfortable handle, are extremely manageable and practical. Choose them according to your interdental spaces; fine for tight spacing, medium for regular spacing and strong for wider spacing.

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
Thanks to its tapered cut, the single tuft brush is used to clean orthodontic braces and their brackets. Its minute head is recommended for plaque removal in interdental spaces, between overlapping teeth and in hard-to-reach areas.

Linear handle designed for easy and correct use of the toothbrush.

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IDEAL FOR ORTHODONTIC BRACES

It can be a real challenge keeping teeth clean when you have braces. With the single tuft brush, it's easier to keep teeth and gums healthy until orthodontic treatment is over, as long as the patient cleans carefully between wires and braces.

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
The extra soft bristles remove food proteins that deposit on the surface of the tongue, neutralising the pathogenic bacteria that cause halitosis.

Starting from the base of your tongue, gently brush with parallel movements. Now brush your teeth with your toothbrush.

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GOODBYE BAD BREATH

In the tongue there are numerous bacteria responsible for bad breath. By eliminating the microflora, therefore, it is possible to reduce bad breath. Cleaning the tongue is also important to slow down the formation of plaque, reduce the risk of gum disease and tooth decay.

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A denture brush allows for effective cleaning of dentures. The dual bristling performs specific dual actions. The large handle allows for controlled brushing, making it easy to clean the denture. Short bristles for cleaning the dental surface. The longer bristles penetrate the hard-to-reach points of the dental arch.

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
The dental mirror allows you to monitor oral health in the hidden areas of the oral cavity, in combination with the bathroom mirror. It is recommended to use it after brushing your teeth well.

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INSTRUCTIONS FOR USE

Thanks to the refractive effect of a second mirror and an adjustable light, it is possible to verify the formation of plaque in the less visible areas of the oral cavity.

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COMPATIBLE ELECTRIC TOOTHBRUSH HEADS

High quality replacement brush heads, compatible with Oral-B* electric toothbrushes. Easy to insert. The rotating head wraps and cleans each tooth and the multiple bristles thoroughly clean the interdental spaces. Dentists recommend changing the head every 2-3 months. Not suitable for children under 4 years.

* Oral-B is a registered trademark of Procter & Gamble.

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SPAZZOLIFICIO PIAVE SpA