

# Galbiati

EXPERIENCE AND QUALITY SINCE 1935

"When a product does not work how it should, or our customers don't like it, for us is a challenge. We feel that from our factory only good products, functional and accessible to all must be produced.

Our daily commitment is to satisfy our customer at 360° and be present before, during and after the purchase of Galbiati products, with the same passion and dedication".



via Concordia 10  
20090 Assago Milano - ITALY

[www.galbiati.com](http://www.galbiati.com)  
[info@galbiati.com](mailto:info@galbiati.com)

## DENTAL SCOUT



# DENTAL UNIT SCOUT

OUTSTANDING DESIGN, GREAT RELIABILITY

## Our Galbiati Lab, 100% Italian

Department of "research and development", created SCOUT, a Premium Class dental unit, entirely designed in Italy, with a great sense of ergonomics and performances.

Solid, easy to use with an uncommon elegance, it's the new born of the GALBIATI family, the first Italian production company of dental units.

## TECHNICAL SPECIFICATION

- **Type of unit:** dental chair with lift system, backrest, motors for chair movements, foot control
- **Supply voltage:** 100 v ac (US) – 230 v ac (EU)
- **Rated frequency:** 50 hz / 60 hz
- **Max power:** 50 hz 650 va
- **Input air pressure:** 0,5 mpa
- **Atmospheric pressure:** 70 kpa to 106 kpa
- **Input water pressure:** 0,6 mpa
- **Water hardness:** 8,4 dh to 12 dh

## PERFORMANCES

- Operating regulation of instruments speed and power
- Reverse of micro motors operating speed
- Adjustment of micro motors speed limit (0-100%)
- Instrument disinfection tubes switch on/off
- Integrated safe system (chair blocked while the instruments are working)
- 4 chair memories
- Instruments cooling switch on/off
- Central suction system
- Water supply system with 2 bottles (disinfection and distilled water)



## Spittoon block

### Suction systems:


- Ejector saliva
- Different suction systems

### Hygiene:

- Clean water bottle system
- Disinfection system
- Both systems

## Button-control panel

- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control (by foot control or manual from keyboard)
- Disinfection program - activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control



## Operating lights

ALYA LED



MAIA LED



Multi functional turning foot control

Monitor

Top mounted system

Low mounted system