



Therapeutic stockings  
segment line

  
**MADE IN ITALY**

## medicaltea® therapeutic line\_\_\_\_\_ p. 4

**COMFORT**  
cotton on the skin

*knee-length, thigh-length, single-leg pantyhose,  
pantyhose, pregnancy pantyhose*

**SOFT**  
microfiber fabric

*knee-length, thigh-length, single-leg pantyhose, pantyhose*

**SIMPLY**  
multifiber fabric

*knee-length, thigh-length, single-leg pantyhose, pantyhose*

**POST-PLASTER  
& POST-TRAUMA**

*knee-length*

## teap® antiembolism line\_\_\_\_\_ p. 10

*Knee-length, thigh-length, single-leg pantyhose*

## medicaltea® orthopedic line\_\_\_\_\_ p. 12

*knee support, ankle support, tight support, elastic elbow,  
elastic back support, post-surgery back support*

## medicaltea® sport line\_\_\_\_\_ p. 16

*knee-length, knee support, ankle support,  
compressive sock*

## medicaltea® travel line\_\_\_\_\_ p. 20

*knee-length*



size table \_\_\_\_\_ p. 22



**TEA** is a company specialized in the production of elastic stockings and segments with quality elastic yarn since 1934.

The variety of products is created in order to satisfy different requirements and pathologies.

Thanks to the ongoing medical research, the technological innovation and the use of certified yarn, it is over 90 years that TEA srl is on market.

**MD**

**MEDICAL DEVICES**





# medicaltea® therapeutic line

Line of healing stockings with decreasing elastic compression for treatment and care of venous stasis.

TEA pays carefully attention to the material of which the stocking is composed of. In fact, the yarn is the essential element: it is selected and certified OEKO TEXT and the supplier companies came from Italy and Europe.

**Medical Tea®** compression stockings are certified and are available in 3 compression classes, according to the pathology:

**KKL1** **KKL2** **KKL3**



### **Class 1** | **KKL 1** 18-21 mmHg

- mild varicose veins
- varicosity with a tendency to mild edema
- varicose veins during and after pregnancy
- functional edema
- teleangectasie
- post-plaster, post-trauma in case of fracture or sprain



### **Class 2** | **KKL 2** 23-32 mmHg

- chronic venous insufficiency
- edema
- superficial phlebitis
- pronounced varicosity
- mild lymphedema
- venous ulcers
- after a surgical or sclerosante treatment

### **Class 3** | **KKL 3**

- lymphedema
- chronic vein thrombosis
- deep vein valve insufficiency
- serious varicose syndrome

# knee-length

Soft upper edge, perfect breathability, reinforced heel, easy to put on with the foot slip provided (only for open toe models).

Available in 4 versions:



**COMFORT**  
cotton  
on the skin

**SOFT**  
microfiber  
fabric

**SIMPLY**  
multifiber  
fabric

**POST-PLASTER  
& POST-TRAUMA**

**Art. 711 KKL 1** (Cotton 21%, Polyamide 53%, Elastomer 26%)

**Art. 722 KKL 2** (Cotton 22%, Polyamide 49%, Elastomer 29%)

**Art. 733 KKL 3** (Cotton 32%, Polyamide 31%, Elastomer 37%)

Toe: open. Color: beige. Length: normal, long.

Plus model: +15% width at the calf.

Made to measure available. Sizes: 1 | 2 | 3 | 4 | 5

**Art. 911 KKL 1** (Polyamide 73%, Elastomer 27%)

**Art. 922 KKL 2** (Polyamide 67%, Elastomer 33%)

Toe: open, closed. Color: beige, black.

Length: normal, long. Sizes: 1 | 2 | 3 | 4 | 5

**Art. 110 KKL 1** (Polyamide 75%, Elastomer 25%)

**Art. 220 KKL 2** (Polyamide 71%, Elastomer 29%)

Toe: open. Color: beige.








Length: normal, long. Sizes: 1 | 2 | 3 | 4 | 5

**Art. 11G KKL 1** (Polyamide 75%, Elastomer 25%)

Toe: open. Color: beige.

Length: normal, long. Sizes: 1 | 2 | 3 | 4 | 5

After fractures or sprains.

COMFORT			SIMPLY
 NORMAL open toe KKL 1 KKL 2 KKL 3	 LONG open toe KKL 1 KKL 2 KKL 3	 PLUS open toe KKL 1 KKL 2	 NORMAL and LONG open toe KKL 1 KKL 2
SOFT			POST-PLASTER & POST-TRAUMA
 NORMAL open and closed toe KKL 1 KKL 2	 LONG open and closed toe KKL 1 KKL 2	 NORMAL open toe KKL 1	

# thigh-length

Elastic bandage with silicone bubble for a better adhesion and breathability.  
Foot slip provided (only for open toe models).



Available in 3 versions:



**COMFORT**  
cotton  
on the skin

**Art. 711 KKL 1** (Cotton 21%, Polyamide 53%, Elastomer 26%)

**Art. 722 KKL 2** (Cotton 22%, Polyamide 49%, Elastomer 29%)

**Art. 733 KKL 3** (Cotton 32%, Polyamide 31%, Elastomer 37%)

Toe: open. Color: beige. Length: normal, long.

Plus model: +15% width at the thigh.

Made to measure available. Sizes: 1 | 2 | 3 | 4

**SOFT**  
microfiber  
fabric

**Art. 911 KKL 1** (Polyamide 73%, Elastomer 27%)

**Art. 922 KKL 2** (Polyamide 67%, Elastomer 33%)

Toe: open, closed. Color: beige, black.

Length: normal, long. Sizes: 1 | 2 | 3 | 4

**SIMPLY**  
multifiber  
fabric

**Art. 110 KKL 1** (Polyamide 75%, Elastomer 25%)

**Art. 220 KKL 2** (Polyamide 71%, Elastomer 29%)

Toe: open. Color: beige.

Length: normal, long. Sizes: 1 | 2 | 3 | 4



COMFORT		
 NORMAL open toe  KKL 1 KKL 2 KKL 3	 LONG open toe  KKL 1 KKL 2 KKL 3	 PLUS open toe  KKL 1 KKL 2
SOFT		SIMPLY
 NORMAL open and closed toe  KKL 1 KKL 2	 LONG open and closed toe  KKL 1 KKL 2	 NORMAL and LONG open toe  KKL 1 KKL 2

# single-leg pantyhose

Adjustable belt with velcro, perfect breathability, reinforced heel, easy to put on with the foot slip provided (only for open toe models).



Available in 3 versions:



**COMFORT**  
cotton  
on the skin

**Art. 711 KKL 1** (Cotton 21%, Polyamide 53%, Elastomer 26%)  
**Art. 722 KKL 2** (Cotton 22%, Polyamide 49%, Elastomer 29%)  
**Art. 733 KKL 3** (Cotton 32%, Polyamide 31%, Elastomer 37%)  
 Toe: open. Color: beige. Length: normal, long.  
 Plus model: +15% width at the thigh.  
 Made to measure available. Sizes: 1 | 2 | 3 | 4

**SOFT**  
microfiber  
fabric

**Art. 911 KKL 1** (Polyamide 73%, Elastomer 27%)  
**Art. 922 KKL 2** (Polyamide 67%, Elastomer 33%)  
 Toe: open, closed. Color: beige, black.  
 Length: normal, long. Sizes: 1 | 2 | 3 | 4

**SIMPLY**  
multifiber  
fabric

**Art. 110 KKL 1** (Polyamide 75%, Elastomer 25%)  
**Art. 220 KKL 2** (Polyamide 71%, Elastomer 29%)  
 Toe: open. Color: beige.  
 Length: normal, long. Sizes: 1 | 2 | 3 | 4



COMFORT		
 NORMAL open toe left/right leg  KKL 1 KKL 2 KKL 3	 LONG open toe left/right leg  KKL 1 KKL 2 KKL 3	 PLUS open toe left/right leg  KKL 1 KKL 2
SOFT		SIMPLY
 NORMAL open and closed toe left/right leg  KKL 1 KKL 2	 LONG open and closed toe left/right leg  KKL 1 KKL 2	 NORMAL and LONG open toe left/right leg  KKL 1 KKL 2



# pantyhose & pregnancy pantyhose

Snug body with gusset, perfect breathability, reinforced heel, easy to put on with the foot slip provided (only for open toe models).



Available in 3 versions:



**COMFORT**  
cotton  
on the skin

**SOFT**  
microfiber  
fabric

**SIMPLY**  
multifiber  
fabric

**Art. 711 KKL 1** (Cotton 21%, Polyamide 53%, Elastomer 26%)  
Toe: open, closed.

**Art. 722 KKL 2** (Cotton 22%, Polyamide 49%, Elastomer 29%)  
Toe: open.

Color: beige. Length: normal, long.

Plus model: +15% width at the thigh.

Made to measure available. Sizes: 1 | 2 | 3 | 4

**Art. 911 KKL 1** (Polyamide 73%, Elastomer 27%)

**Art. 922 KKL 2** (Polyamide 67%, Elastomer 33%)

Toe: open, closed. Color: beige, black.

Length: normal, long. Sizes: 1 | 2 | 3 | 4

**Art. 110 KKL 1** (Polyamide 75%, Elastomer 25%)

**Art. 220 KKL 2** (Polyamide 71%, Elastomer 29%)

Toe: open. Color: beige.

Length: normal, long. Sizes: 1 | 2 | 3 | 4

COMFORT				SIMPLY	
	NORMAL open and closed toe KKL 1		NORMAL open toe KKL 2		NORMAL open toe KKL 1 KKL 2
	LONG open and closed toe KKL 1 KKL 2		LONG open and closed toe KKL 1 KKL 2		LONG open toe KKL 1 KKL 2
	NORMAL open and closed toe KKL 1 KKL 2		LONG open and closed toe KKL 1 KKL 2	<b>COMFORT and SIMPLY</b>	
	PREGNANCY NORMAL open and closed toe KKL 1		PREGNANCY NORMAL open toe KKL 2		



## antiembolism line

Therapeutic stockings are recommended for prevention and care of thromboembolism in lower limbs during and after surgery and in cases of reduced mobility. All the models are manufactured with an inspection hole under toes. The antiembolism line is produced with a particular fabric, LATEX FREE.



# knee-length 18-20 mmHg



**Art. 810 - NORMAL**

**Art. 820 - LONG**

(Polyamide 81%, Elastomer 19%)

*Soft upper edge.*

**Toe:** closed.

**Color:** white.

**Sizes:** S | M | L | XL



# thigh-length 18-20 mmHg

**Art. 811 - NORMAL | Art. 821 - LONG** (Polyamide 81%, Elastomer 19%)

*Elastic bandage with silicone bubble for a better adhesion and breathability.*

**Toe:** closed. **Color:** white. **Sizes:** S | M | L | XL



# single-leg pantyhose 18-20 mmHg

**Art. 812 - NORMAL | Art. 822 - LONG** (Polyamide 81%, Elastomer 19%)

*Adjustable belt with clip.*

**Toe:** closed. **Color:** white. **Sizes:** S | M | L | XL





# medicaltea®


## orthopedic line

This unisex line is particular suitable for the care of knee, ankle, elbow, thigh and lower back.

In this line, TEA uses different types of yarn to create appropriate products which succeeding in fulfilling customers' necessities, like for example: little trauma, edema, joint or post-surgery pain. Articles have to be combined with therapies to be more helpful.

# knee support

UNISEX



## Art. 930 - COTTON

(Cotton 49%, Viscose 20%, Polyamide 16%, Rubber 11%, Elastomer 4%)

*Suitable for little trauma, bruise and post-surgery.*

**Color:** beige. **Sizes:** 1 | 2 | 3 | 4

## Art. 540 - WOOL

(Wool 81%, Cotton 14%, Elastomer 3%, Polyamide 2%)

*Suitable for rheumatic disorders and joint pain.  
It exerts a thermal protective action.*

**Color:** beige. **Sizes:** 1 | 2 | 3 | 4



# ankle support

## Art. 945 - COTTON

(Cotton 47%, Polyamide 36%, Elastomer 17%)

*Suitable for little trauma, edemas,  
sprains, tendinitis.*


**Color:** beige.  
**Sizes:** 1 | 2 | 3 | 4

UNISEX



# thigh support

UNISEX



## Art. 900 - COTTON


(Cotton 52%, Polyamide 28%, Elastomer 20%)

*Suitable for little trauma, support and  
protection.*

**Color:** beige.  
**Sizes:** 1 | 2 | 3 | 4 | 5

# elastic back support

## UNISEX



### Art. 420

(Chlorofiber 61%, Polyamide 16%, Elastomer 12%, Acrylic Fiber 11%)

*It exerts a comfortable thermal protective action.*

**Fabric:** chlorofiber. / **Color:** white.

**Sizes:** 1 | 2 | 3 | 4 | 5



## UNISEX



### Art. 430

(Cotton 73%, Polyamide 15%, Rubber 12%)

*Suitable for post-surgery and post-natal period.*


**Fabric:** cotton, front closure with Velcro.

**Color:** white.

**Sizes:** 1 | 2 | 3 | 4 | 5



## UNISEX



### Art. 440

(Cotton 73%, Polyamide 15%, Rubber 12%)

*Suitable for post-surgery and post-natal period.*

**Fabric:** cotton, front closure with double Velcro.


**Color:** white.

**Sizes:** 1 | 2 | 3 | 4 | 5



# compressive elbow support

UNISEX



## **Art. 901 - SHORT**

(Polyamide 65%, Elastomer 35%)


*A great aid to stimulate the microcirculation.*

**Fabric:** microfiber.

**Color:** beige.

**Sizes:** M | L

UNISEX



## **Art. 902 - LONG**

(Polyamide 65%, Elastomer 35%)

*A great aid to stimulate the microcirculation.*

**Fabric:** microfiber.

**Color:** beige.

**Sizes:** M | L







# medicaltea® sport line

Medical devices that offer a perfect support during physical activity. This line is suitable for all sport types, if used in therapy and prevention.



# knee support

UNISEX

**Art. 300**


(Viscose 62%, Polyamide 20%, Rubber 13%, Elastomer 5%)

*Suitable for support, protection and prophylaxis and treatment of sport rheumatism.*

**Fabric:** viscose.

**Color:** white, blue.

**Sizes:** 1 | 2 | 3 | 4



# ankle support

**Art. 301**

(Cotton 47%, Polyamide 36%, Elastomer 17%)

*Suitable for prevention and prophylaxis and treatment of sport rheumatism.  
Valid support at the ankle.*

**Fabric:** viscose.

**Color:** white, blue.

**Sizes:** 1 | 2 | 3 | 4

UNISEX



## compressive sock KKL.I

### Art. 302

(Polyamide 79%, Elastomer 21%)

*The compression starts already in the plantar arch, it helps to stimulate the microcirculation and to make the sport performance better.*

**Fabric:** microfiber and polyamide.

**Color:** orange/black, yellow/black, green/black, light blue/black.

**Sizes:** 1 | 2 | 3 | 4



## compressive knee-length KKL.I

### Art. 303 (Polyamide 79%, Elastomer 21%)

*The compression starts already in the plantar arch, it helps to stimulate the microcirculation and to make the sport performance better.*

**Fabric:** microfiber and polyamide.

**Color:** orange/black, yellow/black, green/black, light blue/black.


**Sizes:** 1 | 2 | 3 | 4



Made in collaboration with specialists in sport.

# compressive elbow support

UNISEX



## **Art. 901 - SHORT**

(Polyamide 65%, Elastomer 35%)

*A great aid to stimulate the microcirculation.*

**Fabric:** microfiber.

**Color:** black.

**Sizes:** M | L

UNISEX



## **Art. 902 - LONG**

(Polyamide 65%, Elastomer 35%)

*A great aid to stimulate the microcirculation.*

**Fabric:** microfiber.

**Color:** black.

**Sizes:** M | L





## medicaltea® travel line

These compression stockings are useful for those people who attend long flights and have to sit down for hours. These stockings can be the best ally because they reduce the risk of having "traveler's thrombosis", by accelerating blood venous return towards the heart.

# compressive knee-length KKL.I

UNISEX



**Art. 912**

(Polyamide 73%, Elastomer 27%)

*Soft upper edge, perfect wearability.*

**Toe:** closed.

**Length:** normal, long.

**Color:** black, gray.

**Sizes:** 1 | 2 | 3 | 4 | 5




## Size table

medicaTea® therapeutic line  
medicaTea® travel line



### How to identify the size and the length



- Ⓑ ankle (the narrowest point)
- Ⓒ calf (the largest point)
- Ⓓ thigh (under the groin)

#### SIZE

Measure the circumference at the following points:

**B - C** for knee-length

**B - C - G** for thigh-length, single-leg pantyhose, pantyhose, pregnancy pantyhose

*Note:*

*B: circumference of the narrowest point of the ankle*

*C: circumference of the largest point of the calf*


*G: circumference 5 cm under the groin*

#### Circumferences for stockings **NORMAL** model

Size	I	II	III	IV	V
<b>B</b> cm	20-22	23-25	26-28	29-31	32-34
<b>C</b> cm	31-34	36-39	40-44	45-49	50-54
<b>G</b> cm	46-50	52-57	58-63	64-70	

#### Circumferences for stockings **PLUS** model

Size	I	II	III	IV
<b>B</b>	20-22	23-25	26-28	29-31
<b>C</b>	up to cm 37	up to cm 42	up to cm 48	up to cm 53
<b>G</b>	up to cm 55	up to cm 62	up to cm 69	up to cm 77



- Ⓐ heel
- Ⓓ under the knee
- Ⓔ under the groin

#### LENGTH

Measure the length at the following points:

**A - D** for knee-length


**A - G** for thigh-length, single-leg pantyhose, pantyhose, pregnancy pantyhose

#### Leg's length

	NORMAL model	LONG model
<b>A - D</b>	up to cm 39	over cm 39
<b>A - G</b>	up to cm 72	over cm 72



## How to identify the size and the length



- (A) heel
- (B) ankle (the narrowest point)
- (C) calf (the largest point)
- (D) under the knee
- (G) thigh (under the groin)

**SIZE**

Measure the circumference at the following points:

**B - C** for knee-length

**B - C - G** for thigh-length, single-leg pantyhose

*Note:*

*B: circumference of the narrowest point of the ankle*

*C: circumference of the largest point of the calf*

*G: circumference 5 cm under the groin*

## Leg's circumferences

Size	S	M	L	XL
<b>B</b>	19-21 cm	22-24 cm	25-27 cm	28-29 cm
<b>C</b>	30-33 cm	33-38 cm	38-42 cm	42-45 cm
<b>G</b>	up to cm 48	up to cm 55	up to cm 62	up to cm 74

**LENGTH**

Measure the length at the following points:

**A - D** for knee-length


**A - G** for thigh-length, single-leg pantyhose

## Leg's length


	NORMAL model	LONG model
<b>A - D</b>	up to cm 39	over cm 39
<b>A - G</b>	up to cm 72	over cm 72

## How to wear stockings

### OPEN TOE stockings




1 Put on the foot slip.




2 Put on the stockings and remove the slip. Then go ahead as for closed toe stockings.


### CLOSED TOE stockings



1 Roll down the stocking to the heel gently and insert the foot.



2 Pull up the stocking in a regular manner.



3 Continue to pull up the stocking until the position of the stocking is completed.

## ADVISES

- Wear TEA's stocking in the morning, before getting up and take them off before going to bed.
- For a better preservation and duration, TEA recommends a frequent handwash and suggests using products for delicate fabrics. Do not wear the same stockings for over 3 days.
- Take off rings and pieces of jewelry or other stuffs so as to avoid damaging the stocking. Pay attention also at nails, wear gloves instead.
- Check your shoes before putting on the garment and remove anything that can cause damages at the stockings. In case of small ladders, these can be sewed by hand.
- If you notice some threads, do not cut them because they do not represent any mistake.

## CONTRAINDICATIONS

- **Do not use the stocking with damaged skin or deformed limbs.**
- **Do not put in contact with injured skin.**



**Washing instructions:** hand or machine washable with mild detergent.  
Leave to dry away from heat.



## Size table

medicaTea® orthopedic line  
medicaTea® sport line



### KNEE SUPPORT

Measure the circumference of the **knee** to identify the size.

Size	I	II	III	IV
<b>Knee's circumference</b>	30-35 cm	36-39 cm	40-44 cm	45-49 cm



### ANKLE SUPPORT

Measure the circumference of the narrowest point of the **ankle** to identify the size.


Size	I	II	III	IV
<b>Ankle's circumference</b>	19-21 cm	22-24 cm	25-27 cm	28-29 cm



### THIGH SUPPORT

Measure the circumference of the **thigh** to identify the size.

Size	I	II	III	IV	V
<b>Thigh's circumference</b>	42-44 cm	45-50 cm	51-56 cm	57-62 cm	63-70 cm



### ELBOW SUPPORT

Available in **long** or **short** model.

## Size table


medicaTea® orthopedic line  
medicaTea® sport line



### ELASTIC BACK SUPPORT

Measure the circumference of the **waist** to identify the size.

Size	I	II	III	IV	V
<b>Waist's circumference</b>	60-70 cm	70-85 cm	85-100 cm	100-151 cm	151-180 cm



### POST-SURGERY ELASTIC BACK SUPPORT & POST-SURGERY ELASTIC BACK SUPPORT DOUBLE VELCRO

Measure the circumference of the **waist** to identify the size.

Size	I	II	III	IV	V
<b>Waist's circumference</b>	65-80 cm	80-95 cm	95-110 cm	110-125 cm	125-150 cm



### COMPRESSIVE SOCK KKL 1

Measure the circumference of the **ankle** to identify the size.

Size	I	II	III	IV
<b>Ankle's circumference</b>	20-22 cm	23-25 cm	26-28 cm	29-31 cm



### COMPRESSIVE KNEE-LENGTH KKL 1 - SIZE

To identify the size, measure at the following points:

Size	I	II	III	IV
<b>ANKLE</b> (the narrowest point)	20-22 cm	23-25 cm	26-28 cm	29-31 cm
<b>CALF</b> (the largest point)	31-34 cm	36-39 cm	40-44 cm	45-49 cm

### COMPRESSIVE KNEE-LENGTH KKL 1 - LENGTH

To identify the length, measure at the following points:

Size	NORMAL model	LONG model
<b>from HEEL to UNDER KNEE</b>	up to 39 cm	over 39 cm



It is over 90 years that TEA has passion and experience in the field of therapeutic stockings and segment line.

All the manufacturing cycle is entirely produced in TEA's textile plant in Avigliana (TO - Italy).

It starts from the spools of yarn, then there is the design of products, the production, the sewing, the ironing, the measurement of sizes and compression, the quality control, the labelling and in the end there is the packaging, where the product is ready to be sent to customers.

TEA Quality System ISO 9001 and ISO 13485 is certified by Bureau Veritas.





**TEA S.r.l.**

Via A. Einstein, 17 - 10051 Avigliana (TO) Italy  
Tel. +39 011 9348971 - [medical@medicaltea.it](mailto:medical@medicaltea.it)  
[www.medicaltea.it](http://www.medicaltea.it) - [www.calzemonouso.it](http://www.calzemonouso.it)  
[www.medicalteashop.it](http://www.medicalteashop.it)



**MD** Medical Device Reg.(UE)  
2017/745 and s.m.s.

TEA Quality System ISO 9001  
and ISO 13485 is certified  
by Bureau Veritas.



   
**MADE IN ITALY**

**ISTRUZIONI PER LO SMALTIMENTO**

Il prodotto deve essere smaltito  
in conformità alle normative vigenti



**OPUSCOLO  
PAP 22**

**► RACCOLTA  
CARTA**

Verifica le disposizioni  
del tuo Comune.

